EYE MAX 2 MICRO

EN MICRO: OFC 20 - 25 CM (8 - 10 IN)

ENGLISH - Indications for Use:

The EyeMax2 provides neonatal eye protection during ultraviolet (UV) phototherapy treatment of jaundice. The EyeMax2 conforms with EN 170:2002 Clause 5.2 (Scale Number 2-5) for spectral transmittance as appropriate for the indication for use.

Instructions for use:

- Measure Occipital-Frontal Circumference (OFC) of the patient.
- Select the appropriate EyeMax2 size based on the OFC. See chart for appropriate sizing.
- Disconnect the Velcro.
- Position the eye pads over the eyes.
- Pass the lower strap under the small of the neck and attach velcro.
- Secure the Velcro of both straps and position the EyeMax2 as shown.
- Make sure the strap completely covers the ears.
- If the eye patch rides up away from the nose, tighten the lower strap.
- If the eye patch rides down towards the nose and mouth, repeat application process beginning with step 1.

SYMBOL GUIDE

The following symbols and safety labels are found on the EyeMax2:

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
EyeMax2, strana 1

EyeMax2, strana 2

MD	Medical Device	Rxonly	RX Only
UDI	Unique Device Identifier	UK RP	Responsible Person in the UK

⚠️ Precautions:

CAUTION: Improper EyeMax2 sizing may lead to patient airway obstruction or inadequate eye protection. The Eye max2 should only be used while the patient is under the continuous, direct supervision of healthcare professionals.

- ⓧ DO NOT** tug on the EyeMax2 to reposition it. A premature infant's skin is very delicate and can easily rupture. Avoid obstructing the nose. Check regularly to ensure the EyeMax2 is not blocking the patient's airway. The EyeMax2 is not a bandage. Cover any existing skin abrasions with appropriate bandages before use.
- ⓧ DO NOT** over tighten - causing excessive pressure on the eyes. Make sure the strap is covering the ear, not hooked over the ear, and that the ear is not folded over. Replace the EyeMax2 with a new one if it becomes over-stretched. Reuse of this device could cause cross-contamination or product malfunctions.

R300P01 Regular: occipital-frontal circumference 33-38 cm in packs of 20.

R300P02 Preemie: occipital-frontal circumference 26-32 cm in packs of 20.

R300P03 Micro: occipital-frontal circumference 20-25 cm in packs of 20.

ⓧ	Single Patient Use Only Non-Sterile, Latex-Free
----------	--

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EyeMax2, strana 69

EyeMax2, strana 70

EN: Note: Serious incident(s) that occur in relation to the device should be reported to Maxtec and the competent authority of the Member State in which the user and/or patient is established. Serious Incident(s) is defined as directly or indirectly led, might have led, or might lead to the death of a patient, user, or other person; the temporary or permanent serious deterioration of a patient's, user's or other person's state of health; or a serious public health threat.

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