KONICA MINOLTA

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1. DESCRIPTION OF THE ANALYZER / PARTS

For the following descriptions, please refer to the diagrams on the inside back cover of this instruction manual. See **A1** for names of parts, **A2** for Display, and **A3** for Probes.

- A1**
- A 1 Power Button (to switch ON/OFF)
 - 2 LCD (to display the measurement value and messages)
 - 3 Backlight button
 - B 4 Probe connector (to connect the probe)
 - 5 External output terminal
 - 6 Battery cover (to lock the chamber)
 - 7 Wrist band (to attach the instrument on the wrist)
 - 8 Band fixture A
 - 9 Band fixture B
- A2**
- 1 Pulse rate (P.R.)
 - 2 SpO₂ value
 - 3 Battery mark (indicates the remaining battery level in 3 steps)
 - 4 Pulse level Meter
- A3**
- 1 Finger Clip Probe SR-5C
 - 2 Spot Check Probe SP-5C
 - 3 Monitor Probe LM-5C
 - 4 Personal Probe SD-5C
 - 5 Extension Cable EC-300

Calculation Method for Displayed Values

<Pulse level meter>

The pulse level meter displays the transmitted light amount ratio calculated by the following equation. The pulse level is indicated in eight levels (0 to 15%, full scale if over 10%).

$$\text{Pulse level(\%)} = \frac{\text{Transmitted light amount that varies due to the pulse}}{\text{Transmitted light amount that does not vary due to the pulse}} \times 100$$

<Oxygen saturation (SpO₂)>

Updates the moving average of the oxygen saturation over the last 12 seconds at one-second intervals.

<Pulse rate (P.R.)>

Updates the moving average of the last 8 pulse rates at one-second intervals.

2. PREPARATIONS

Before starting measurements using this instrument, follow the procedure given below to ensure that measurements are stable.

1) Attaching the Wrist Band

Attach the Wrist Band WB-300 (supplied as a standard accessory) to the instrument. Since the WB-300 is made of elastic material, it allows for flexibility and comfort of fitting on the wrist.

Procedure

- 1 Place the instrument with its rear side up.
- 2 Pass the angled tip of the hook-and-loop tape of the wrist band outward through band fixture A.
 - ◆When doing so, make sure that the Velcro tape section faces upward.
- 3 Fold up the inserted hook-and-loop section, and stick it firmly to the band.
- 4 Pass the other end of the band outward through band fixture B.
 - ◆Take care not to twist the band.
- 5 Fold back and secure the band with the hook-and-loop tape section.
 - ◆When attaching the instrument to an arm, adjust the band length and fitting at this hook-and-loop tape section.

2) Connecting the probe

Connect the probe's connector plug into the Probe Connector on the instrument. When using the Extension Cable EC-300, connect the extension cable EC-300 between the PULSOX-300i and the probe.

It extends the probe cable to a total length of 1m.

- ◆ Make sure that probe's connector plug is oriented correctly.
- ◆ Do not connect it in with excessive force. Doing so may damage the connector.
- ◆ Do not touch the connector terminals by hand. Contact failure may occur.

3) Installing the battery

Procedure

- 1 Turn the battery cover clockwise approximately 90 degrees to open the cover.
 - ◆ When turning the cover, keep the flat part of your finger on the cover's catch. This facilitates opening the cover with minimum force.
- 2 Insert the battery into the battery compartment according to the polarity mark(⊕) shown on the rear of the instrument.
 - ◆ This instrument requires one AAA-size battery.
- 3 Turn the battery cover counter-clockwise approximately 90 degrees to close it.
 - ◆ Introduction of a new mechanism locks the cover automatically when it is closed.

3. MEASUREMENT

1) Starting the Measurement

- 1 Attach the probe to a finger of an adult patient.
- 2 Press the Power button for a long time.
 - The Power is switched ON, all the elements on the display are displayed for about 2 seconds, the set data in the instruments (day and time, next measurement file number, etc.) are displayed for each 2 seconds, and about 10 seconds after the power is switched ON, measurements will be started.
 - ◆ Until the measurement value has been calculated and displayed, "---" will be displayed as the values of the oxygen saturation (SpO₂) and the pulse rate.

2) Reading the display

When the instrument is turned ON, all the elements of the LCD will light up. The pulse level, pulse rate and SpO₂ value are displayed during measurement.

3) Ending the Measurement

Press the Power button for a long time, the measurement is finished and the Power is off.

About Auto Power-Off function

By connecting the instrument to a computer and using the optional software, the time to be elapsed from power-on (start of measurement) to power-off (end of measurement) can be set as the auto power-off timer. If an auto power-off timer is set, the measurement is finished automatically in the time after power-on.

Regarding the memory function

- The memory built into the PULSOX-300i will continue to hold data even if the battery is removed.
- A maximum of 300 hours of measurement data can be stored; i.e., 10 hours of measurement data can be stored for 30 persons. Storage for measurement data is set up in 30-minute blocks, so even if the duration of measurement data is less than 30 minutes, it will be treated as a 30-minute block. For this reason, if measurement data with durations of shorter than 30 minutes are stored, the actual measurement duration will be different from the stored duration, and the maximum storable time will be shorter than 300 hours.
- The maximum storable number of data files is 399 files. After that, the No. 400 data file will be stored as data file No. 1, overwriting the previous No. 1 data file.



4. NOTES ON USE




- This instrument should be used under the following operating conditions.
 - Temperature of 0 to 40°C (32 to 104°F), relative humidity of 30 to 85%, with no condensation.
 - Atmospheric pressure of 700 to 1060 hPa (altitude of -400 to 3000m).




- This is a precision instrument. To avoid the possibility of it being damaged, the instrument should not be dropped nor should heavy objects be placed on top of it.
- Do not expose this instrument to rain and water.
- This instrument has no alarm function. Do not use this instrument if an alarm function is necessary.
- This instrument is designed for use on adults. When using the instrument on infants or babies, please note that the designed measurement accuracy cannot be guaranteed.
- Do not use the instrument on a point where there is excessive vibration. In addition, to avoid excessive impact on the instrument, it should be handled gently. Failure to observe this may cause breakdown.
- A maximum of 300 hours of measurement data can be stored, or the maximum storable number of data files is 399 files. If measurements are continued after the maximum number of storable data has been reached, the new data will overwrite the oldest data. To keep data, be sure to download it before the memory becomes full.
- Accurate measurement data may not be obtained in the following cases.
 - When body movement is excessive
 - When the probe is not attached properly
 - When blood circulation at the measuring point is poor (due to pressure on arm/finger, peripheral circulatory insufficiency)
 - When ambient light (e.g. panel light, fluorescent light, infrared heat lamp, direct sunlight) is too strong
 - When influenced electromagnetically by other electronic devices (e.g. near electrical appliances like TV, medial device)
 - When a mobile telephone is used during measurement
 - When influenced by abnormal hemoglobin like carbon monoxide hemoglobin (HbCO) and methemoglobin
 - When pigments such as cardiogreen, intravascular dyes and indocyanine green are present in the blood
 - When finger nails are polished




5. ERROR MESSAGES


The instrument displays error messages when the error occurs. Check and correspond to them correctly.

Error Messages	Cause	Solution
 <p>• "C" blinks.</p>	<ul style="list-style-type: none">• There is no probe connected to the PULSOX-300i or the probe is improperly connected to the PULSOX-300i.• If using the optional Extension Cable EC-300, the cable may not be connected to the instrument or the probe correctly.	<ul style="list-style-type: none">• Connect the probe to the PULSOX-300i properly. When using the optional Extension Cable EC-300, connect the cable to the instrument and the probe correctly.• If this message reappears even though the probe is properly connected to the PULSOX-300i, the probe is malfunctioning (e.g., wire breakage inside the probe). Use another probe.
 <p>• "L" blinks.</p>	<ul style="list-style-type: none">• Insufficient light for measurement.• The light emitted from the LED enters the sensor directly.• Strong light enters the probe's sensor directory.	<ul style="list-style-type: none">• Check that the probe is properly attached to the patient.• Be sure that the light-source and sensor sections are properly lined up.• If the part being measured is thick, insufficient light will be transmitted. Attach probe to a thinner part.• Clean the light-source and sensor sections.• Take measures to prevent direct entry of strong light to the probe.• If this message reappears even though the above points have been checked, the probe may be malfunctioning. Use another probe.

Error Messages	Cause	Solution
 <ul style="list-style-type: none"> • Measurement value and "L" flash alternately. 	<ul style="list-style-type: none"> • The light power required for the measurement are not obtained. The measurement value is displayed, however, the accuracy of the measurement values cannot be guaranteed. 	<ul style="list-style-type: none"> • Check that the probe is properly attached to the patient. • Be sure that the light-source and sensor sections are properly lined up. • If the part being measured is thick, insufficient light will be transmitted. Attach probe to a thinner part. • Clean the light-source and sensor sections. • Take measures to prevent direct entry of strong light to the probe. • If this message reappears even though the above points have been checked, the probe may be malfunctioning. Use another probe.
 <ul style="list-style-type: none"> • "P" blinks. 	<ul style="list-style-type: none"> • Pulse signal required for measurement is not being received. 	<ul style="list-style-type: none"> • Check that the probe is properly attached to the patient. If the probe is properly attached but this message still appears, the circulation in the area being measured may not be good. Move the probe to an area from which a sufficiently strong pulse signal can be received. • If a sufficiently strong pulse signal cannot be obtained regardless of the area to which the probe is attached, the probe may be malfunctioning. Use another probe.
 <ul style="list-style-type: none"> • Measurement value and "P" flash alternately. 	<ul style="list-style-type: none"> • The pulse is weak and the pulse signals required for the measurement are not obtained. The measurement value is displayed, however, the accuracy of the measurement values cannot be guaranteed. 	<ul style="list-style-type: none"> • Attach the probe to the body correctly. If this message is displayed when the unit is attached to the patient correctly, the measurement conditions are regarded as blood circulation or other circulatory problems. If this occurs, search for a location where adequate pulse signals can be obtained. Attach the probe again or warm up the measurement location. • In particular, when the pulse is weak, this error tends to occur. For improved measurement accuracy, if the probe is exposed to strong light, wrap it with a piece of black cloth, etc. Also make sure that the patient rests calmly to avoid influence of body movement.

Error Messages	Cause	Solution
 <ul style="list-style-type: none"> • "H" blinks. 	<ul style="list-style-type: none"> • The pulse is too strong for correct measurement. 	<ul style="list-style-type: none"> • Attach the probe to the body correctly. If this message is displayed when the unit is attached to the patient correctly, the measurement conditions are regarded as blood circulation or other circulatory problems. If this occurs, search for a location where more stable pulse signals can be obtained. • If this message reappears even though the above points have been checked, the probe may be malfunctioning. Use another probe.
 <ul style="list-style-type: none"> • Measurement value and "A" flash alternately. 	<ul style="list-style-type: none"> • Motion artifact The measurement value is displayed, however, the accuracy of the measurement values cannot be guaranteed. 	<ul style="list-style-type: none"> • Attach the probe to the patient properly. If this message still reappears even though the probe is attached properly, keep the measuring point as stationary as possible or attach the probe to another part of the body where body movement is relatively low.
 <ul style="list-style-type: none"> • "O-" blinks. 	<ul style="list-style-type: none"> • The pulse rate exceeds the upper limit (230 bpm) of the measurement range. 	<ul style="list-style-type: none"> • Measurement cannot be performed when the upper limit of the measurement range is exceeded.

Error Messages	Cause	Solution
 <ul style="list-style-type: none"> • "-U-" blinks. 	<ul style="list-style-type: none"> • The pulse rate belows the lower limit (30 bpm) of the measurement range. 	<ul style="list-style-type: none"> • Measurement cannot be performed when the lower limit of the measurement range is belowed.
 <ul style="list-style-type: none"> • Measurement value blinks. 	<ul style="list-style-type: none"> • The oxygen saturation (% SpO₂) is alue below the lower limit (70%) of the measurement range. The measured value will still be displayed, but its accuracy cannot be guaranteed. 	<p>(This indicates that the measurement value is out of the guaranteed accuracy range.)</p>
	<ul style="list-style-type: none"> • Battery Out 	<ul style="list-style-type: none"> • Replace the battery with new one. • In this case, the measurement data is not left if the battery is removed.
 <ul style="list-style-type: none"> • Battery mark is displayed as level 2. 	<ul style="list-style-type: none"> • Battery power is low. 	<ul style="list-style-type: none"> • For about 8 hour's measurement is available. • If it is necessary to take measurement over 8 hours, remove battery and replace with new one. • If the battery is out in the middle of measurement, the measurement value is not left.
 <ul style="list-style-type: none"> • Battery mark is displayed as level 3 only. 	<ul style="list-style-type: none"> • Battery power is empty. 	<ul style="list-style-type: none"> • Replace the battery with new one as soon as possible. • If the battery is out in the middle of measurement, the measurement value is not left.

Error Messages	Cause	Solution
	<ul style="list-style-type: none"> Any trouble may occur in the measurement circuit of the instrument. 	<ul style="list-style-type: none"> Turn the Power Off, and turn ON again. If these message reappear even though the above points have been checked, the instrument may be malfunctioning. Contact to the nearest distributor or authorized service facilities.

6. SAFETY PRECAUTIONS

To ensure correct use of this instrument, read the following points carefully and adhere to them. After you have read this manual, keep it in a safe place where it can be referred to anytime a question arises.

Warning

(Failure to adhere to the following points may result in death or serious injury.)

- The instrument is designed for measurement of the oxygen saturation (SpO₂) of arterial blood and the pulse rate. Do not use it for any other purposes, such as warning of sleep apnea and breathing abnormalities.
- Do not use the instrument in places where flammable or combustible gases (anesthetic gas etc.) are present. Doing so may cause a fire.
- Do not disassemble or modify the instrument and accessories. Doing so may cause a fire or electric shock.
- The instrument should not be operated if it is damaged, or smoke or odd smells occur. Doing so may result in a fire. In such situations, turn OFF the power immediately, remove the battery, and contact the nearest authorized service facility.
- Do not put the batteries on a fire, short-circuit them, heat them or disassemble them. Doing so may cause explosion or heat generation, resulting in fire or injury.

Caution

(Failure to adhere to the following points may result in injury or damage to the instrument or other property.)

- Do not use batteries other than those specified by KONICA MINOLTA SENSING, INC. When installing batteries in the instrument, make sure that they are correctly oriented according to the mark.
- If alkali fluid from the battery comes in contact with eyes, skin, or clothing, immediately wash the affected area and see a physician for treatment.
- Do not use probes other than those specified by KONICA MINOLTA SENSING, INC. Use of other probes may cause the probe to overheat, resulting in burns.
- Do not operate the instrument for long periods of time with a probe

attached to a patient. Low-temperature burn, redness or rash may result. If you feel pain or itchiness, stop use of the instrument immediately and consult a doctor. A doctor should also be consulted before using the instrument on patients with high fever, those with peripheral blood circulation problems or those with sensitive skins.

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© 2006 KONICA MINOLTA SENSING, INC.