

MaxVenturi Usage Sheet

SENSING

ANALYSIS

DELIVERY

Approximate MaxVenturi Oxygen Consumption (LPM)

		<i>Flow Setpoint</i>					
		10	20	30	40	50	60
O₂%	32%	1.4	2.8	4.2	5.6	7.0	8.3
	40%	2.4	4.8	7.2	9.6	12.0	14.4
	50%	3.7	7.3	11.0	14.7	18.4	22.0
	60%	4.9	9.9	14.8	19.7	24.7	29.6
	70%	6.2	12.4	18.6	24.8	31.0	37.2
	80%	7.5	14.9	22.4	29.9	37.3	44.8
	90%	8.7	17.5	26.2	34.9	43.7	52.4
	100%	10.0	20.0	30.0	40.0	50.0	60.0

Approximate MaxVenturi Run Time (Hours)

		<i>Flow Setpoint</i>					
		10	20	30	40	50	60
O₂%	32%	7.1	3.6	2.4	1.8	1.4	1.2
	40%	4.2	2.1	1.4	1.0	0.8	0.7
	50%	2.7	1.4	0.9	0.7	0.5	0.5
	60%	2.0	1.0	0.7	0.5	0.4	0.3
	70%	1.6	0.8	0.5	0.4	0.3	0.3
	80%	1.3	0.7	0.4	0.3	0.3	0.2
	90%	1.1	0.6	0.4	0.3	0.2	0.2
	100%	1.0	0.5	0.3	0.3	0.2	0.2

MAXVENTURI USAGE CALCULATOR

The graphs above can help to plan and understand the oxygen consumption when using a MaxVenturi with oxygen coming from an E-cylinder.

Note: *These values are intended as a general guideline only and are not to be interpreted as specifications. Testing in your clinical application is required. Times are based on an estimated 600 Liter 'E' cylinder.*