Blender Buddy

BIRD® BLENDER

EN	Instructions for Use 2	NL	Gebruiksaanwijzingen 8
FR	Instructions d'utilisation . . 4	IT	Istruzioni per l'uso 10
DE	Bedienungsanleitung 6	ES	Instrucciones de uso 12

NOTE: The latest edition of this operating manual can be downloaded from our website at www.maxtec.com

WARRANTY

2-year limited warranty on manufacturing defects.

WARNINGS

-  **DO NOT** use substitute parts.
-  **DO NOT** use lubricants on o-rings.
- ◆ Not for use in MRI environments.
- ◆ Carefully read instructions provided with flow meters, as well as this insert, before use.
- ◆ For use with air/oxygen blenders only.
- ◆ Max pressure is 100 psi.
- ◆ Check for leaks and proper function before using on a patient.
- ◆ Maxtec assumes no responsibility for any damage or injury caused by improper installation, assembly or use of this product.
- ◆ Replace EPDM o-rings during recommended blender rebuild. Use Maxtec approved o-rings.
- ◆ This product should only be used under proper supervision of a healthcare professional.
- ◆ If components are damaged or missing, contact your dealer immediately.
- ◆ Clean often using a cloth moistened with 65% isopropyl alcohol/water solution or germicidal wipe.
- ◆ Follow all recommendations provided in the manual for the connected blender.

1.0 INTRODUCTION

The Blender Buddy is a manifold accessory designed to be attached to an outlet port of the Bird® air/oxygen blender. From the blender outlet port, the air/oxygen mixture passes into a primary or a primary and a secondary channel within the Blender Buddy manifold. The primary channel is used to deliver the gas to a set of flowmeters which meter the gas flow by means control valves. The secondary channel is a small channel that empties into a cavity (port) sized to receive a galvanic oxygen sensor; this port fits the Maxtec line of oxygen analyzers and is used to analyze oxygen concentration of the air/oxygen mixture to be delivered to a patient. The Blender Buddy can be configured with various sized color coded flowmeters for customized flow rate delivery to patients. Flowmeters are sold separately.

1.1 Indications For Use


The Blender Buddy is intended for use in conjunction with a Bird® air/oxygen blender. Then Blender Buddy allows the operator of a blender to supply a mixed gas through a set of flowmeters. Blender Buddies designed with bleeds include an extra analyzer port for analyzing the oxygen concentration of the air/oxygen mixture being delivered to a patient.

1.2 Symbol Guide



	Federal law (USA) restricts this device to sale on or by the order of a physician		Follow instructions for use.
	Authorized Representative in the European Community		Do Not
	Manufacturer		Contains no Polyvinyl Chloride
	Date of Manufacture		Medical Device
	Storage Temperature Range		Not Manufactured with Natural Rubber Latex
	MR Unsafe		Warning
	Catalog Number		Lot Number
	Made in USA		

2.0 INSTRUCTIONS FOR USE

Before assembly verify that the flow meters are the correct flow and gas type according to the labeling on the flow meters. Prior to assembling the Blender Buddy inspect all parts for signs of damage.

1. Be sure the small o-ring  is present around the threads of the short bolt  between the small block  and the blender .
2. Attach the small adapter block  to the left side of the blender using the short bolt  and provided Allen Wrench (No thread tape is necessary). Remove any fittings from the left side outlet if necessary.
*Skip to Step 9 if flow meters are already attached.
3. Wrap both threaded ends of the pipe nipples  with thread tape (1-1/2 wraps). Use caution to ensure tape does not occlude the fittings.
4. Use a 7/16" wrench to tighten the pipe nipples  into the threaded holes of the flow meters   .
5. Tighten one assembled flow meter  into one of the threaded holes on the large block  . Use caution to ensure the flow meter is mounted vertically.
6. Remove the clear flow tube and shield from the other flow meter  . Note how the flow meter components are arranged for reassembly in a subsequent step. Use caution to prevent contamination from entering the disassembled flow meter.
7. Tighten the disassembled flow meter  into the other threaded hole in the large block  .
8. Reassemble the flow meter  . Use caution to ensure the flowmeter is mounted vertically.
9. Insert large o-ring  into o-ring groove of the large block  .

10. Verify white bleed control plug is inside port of the large block.
11. Attach large block (A) to small block (B) with the two long screws (C).
12. Supply gas to the blender and check for leaks and proper function.

NOTE: The blender buddy is designed for use on the primary port (left side) of a blender. The blender buddy has a 3 lpm continuous bleed built into the manifold of the device. When properly installed on the left side of a blender the blender buddy's bleed ensures accuracy of flow from the primary port at low flow rates, i.e. less than 3 lpm and is equal to the low flow accuracy of the auxiliary port on the right side of the blender.

For use with the following models:

- Bird® 3800 Microblender
- Bird® Low Microblender
- Bird® Omniblender

2.1 Flow meter options

0 - 200 ccm Flowmeter (white).....	RP34P03-103
0 - 1 LPM Flowmeter (white).....	RP34P03-104
0 - 8 LPM Flowmeter (white).....	RP34P03-105
0 - 15 LPM Flowmeter (white).....	RP34P03-101
0 - 70 LPM Flowmeter (white).....	RP34P03-106
0 - 3.5 LPM Flowmeter (white).....	RP34P03-102
2 - 26 LPM Flowmeter (white).....	RP34P03-108
0 - 5 LPM LPM Flowmeter (white).....	RP34P03-109
0 - 200 ccm Flowmeter (green).....	RP34P03-003
0 - 1 LPM Flowmeter (green).....	RP34P03-004
0 - 8 LPM Flowmeter (green).....	RP34P03-005
0 - 15 LPM Flowmeter (green).....	RP34P03-001
0 - 70 LPM Flowmeter (green).....	RP34P03-006
0 - 3.5 LPM Flowmeter (green).....	RP34P03-002
2 - 26 LPM Flowmeter (green).....	RP34P03-008
0 - 5 LPM Flowmeter (green).....	RP34P03-009

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