Regular: OFC 33 - 38 cm (13 - 15 in)

Indications for Use:

The Eyemax2 provides neonatal eye protection during ultraviolet (UV) phototherapy treatment of jaundice. The Eyemax2 conforms with EN 170:2002 Clause 5.2 (Scale Number 2-5) for spectral transmittance as appropriate for the indication for use.

Instructions for use:

1. Disconnect the Velcro.
2. Position the eye pads over the eyes.
3. Pass the lower strap under the small of the neck and attach velcro.
4. Secure the Velcro of both straps and position the EyeMax2 as shown.
5. Make sure the strap completely covers the ears.
6. If the eye patch rides up away from the nose, tighten the lower strap.

Catalog # R300P01 **Regular:** occipital-frontal circumference **33-38** cm in packs of 20.

Catalog # R300P02 **Premie:** occipital-frontal circumference **26-32** cm in packs of 20.

Catalog # R300P03 **Micro:** occipital-frontal circumference **20-25** cm in packs of 20.

⚠️ Precautions:

- ❌ **DO NOT** tug on the EyeMax2 to reposition it.
A premature infant's skin is very delicate and can easily rupture.
Avoid obstructing the nose. Check regularly to ensure the EyeMax2 is not blocking the patient's airway. The EyeMax2 is not a bandage. Cover any existing skin abrasions with appropriate bandages before use.
- ❌ **DO NOT** over tighten - causing excessive pressure on the eyes.
Make sure the strap is covering the ear, not hooked over the ear, and that the ear is not folded over.
Replace the EyeMax2 with a new one if it becomes over-stretched. Reuse of this device could cause cross-contamination or product malfunctions.

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SGS United Kingdom Limited