

Blender Maintenance Kit



Maxtec's Blender Maintenance Kit is the only comprehensive kit with known compatibility to service the Bird® MicroBlender, Sentry 2 Blender, and most MaxBlend™ 2 Blenders.*


This kit provides all the commonly serviced components required for biennial preventative maintenance for both high and low flow versions of new and existing blenders in the market.

BLENDER MAINTENANCE KIT PART NUMBER

Blender Maintenance Kit (International) 10003V2

**Please refer to blender service manual to ensure kit compatibility*

KIT CONTENTS	QTY
O-RING, 0.117 ID X 0.040	2
O-RING, 0.176 ID X 0.070	10
O-RING -O10, .239 ID X .070	1
O-RING, EPDM 0.364 ID X .504 OD	3
O-RING, 0.114 ID X 0.070	1
O-RING, 0.301 ID X 0.070	1
O-RING, 0.426 ID X 0.070	2
RING, INTERNAL RETAINING 0.39	1
MUFFLER, BLENDER FOAM, SINGLE	1
O-RING, .468 ID X .078	4
ASSEMBLY, BALANCE BLOCK DIAPHRAGM	2
VALVE, DUCKBILL CHECK	2
WASHER, 0.5 OD STEP	2
PLUG, ALARM CAP FOAM	1
O-RING, 0.426 ID X 0.040	2
BALL, CHECK-RUBBER, 0.218 Ø	1
SPRING, 0.093 OD X 0.053 ID X 0.300	1
O-RING, 0.414 ID X 0.070	7
O-RING, 0.614 ID X 0.070	1
O-RING, 0.239 ID X 0.070	3
FILTER, NYLON CONE	2
O-RING, 0.313 ID X 0.051	1
O-RING, 0.410 ID X 0.050	2
O-RING-O12, 0.364 ID X .070	1



Discover more about the
Blender Maintenance Kit

