

Provita® I.V. Poles


SENSING

ANALYSIS

DELIVERY



Regular I.V. Pole



Pediatric I.V. Pole



Folding I.V. Pole
With Bag

I.V. POLE SOLUTIONS

The Provita® Regular I.V. Poles from Maxtec allow a low-cost alternative for patient care stands.

The Pediatric pole combines a colorful presence with a creative design to help create a more comfortable setting in the patient-area.

The Folding stand is ideal for anywhere space-savings during storage is a concern, or for Homecare applications.

PART NUMBERS

Normal I.V. Pole R100P85-001
Pediatric I.V. Pole R100P85-002
Folding I.V. Pole R100P85-003

Normal, Specifications

| | |
|-----------------------------|-----------------------------------------------------------------------------------------|
| Material: | Chromed steel |
| Color of plastic parts: | Blue |
| Bottle holder: | Plastic 4 safety bottle hooks - 2 kg / hook |
| Adjustment: | Screw safety height adjustment |
| Regulating range: | 1.350 - 2.150 mm |
| Tube diameter: | 25 / 18 mm |
| Length of outer tube: | 1.000 mm |
| Base diameter: | 635 mm |
| Base leg profile: | 25 x 25 mm |
| Casters: | Twin casters Ø 50 mm, according to DIN EN 60601 electrically conductive and lockable |
| Total- / base weight (net): | Appr. 3,5 kg / 1,9 kg |



Pediatric, Specifications

| | |
|----------------------------|---------------------------------------------|
| Material: | Aluminum / plastic |
| Color of plastic parts: | Yellow / blue / green |
| Bottle holder: | Plastic 4 safety bottle hooks - 2 kg / hook |
| Adjustment: | Screw safety height adjustment |
| Regulating range: | 1.350 - 2.000 mm |
| Tube diameter: | 25 / 18 mm |
| Base diameter: | 635 mm |
| Casters: | Twin casters Ø 50 mm with airplane casters |
| Total / base weight (net): | Appr. 3,0 / 1,3 kg |



Folding, Specifications

| | |
|----------------------|-------------------------------------------------|
| Material: | Aluminum / plastic |
| Bottle holder: | 2 safety bottle hooks - 0,5 kg / hook |
| Adjustment: | Screw safety height adjustment |
| Regulating range: | 830 - 1.780 mm |
| Tube diameter: | 20 / 15 / 11 mm |
| Base diameter: | 680 mm |
| Total load capacity: | 3 kg |
| Total weight (net): | appr. 2 kg |
| Features: | Anti-slide stopper base, including carrying bag |



ML# 280 REV. A