


Maxtec Serial Number Chart

SENSING

ANALYSIS

DELIVERY



► Maxtec Serial Number Definition

Maxtec Oxygen Sensor serial numbers are located on the oxygen sensor label. The Lot Number, Job Number, Line Item Number and Sequence Number all relate to the sequence of when the sensor was manufactured. The model number is listed separately on the label.

To Determine Warranty Period

Look at the product's serial number. The first two letters indicate the year and month in which the lot was started. (e.g., the date code referenced in the serial number above (ZB) was started in 2014, during the month of February). Warranty periods for each product can be found on our website www.maxtec.com.

ML-0020 REV M