

SmartStack™ Comparison Chart

24" BASE	24" BASE	24" BASE	24" BASE
48 Month Warranty	NA Month Warranty	NA Month Warranty	NA Month Warranty
SMARTSTACK BY MAXTEC	MEDLINE MPH08085.86	WT FARLEY IVS-5110	IV LEAGUE 1085-2/1085-4
Five 6" (two locking) Casters	Five 3" Casters	Five or Six 3" (two locking) Casters	Five or Six 3" Casters
65" to 88" Pole Length	52" to 94" Pole Height	NA Pole Length	55" to 97" Pole Length
Stainless Steele Pole Material	Stainless Steele Pole Material	Stainless Steele Pole Material	Chrome Plated Pole Material
Poweder Coated Aluminum Base Material	Cast Aluminum, Epoxy Coated Base Material	NA Base Material	Powder Coated Base Material
Easy, One-handed Height Adjustment	Manual twist-screw Height Adjustment	Manual twist-screw Height Adjustment	Squeeze Release Height Adjustment
Yes Stackable	No Stackable	No Stackable	No Stackable

*Maxtec is not affiliated with nor endorsed by Medline, WT Farley or IV League.



2305 South 1070 West
Salt Lake City, Utah 84119

866.4.Maxtec | www.maxtec.com