

December 11th, 2019

Rail Mount Bracket Discontinuation

Dear Valued Customer,

This letter is to inform you of a change to a product(s) you've previously purchased from Maxtec. We are removing the blender and humidifier rail mount brackets from our system, as we now offer a more convenient bracket that is compatible to mount to both poles and rails.

PART NUMBERS BEING REMOVED

- BLENDER RAIL MOUNT BRACKET **Part Number: R100P09**
- HUMIDIFIER RAIL MOUNT BRACKET **Part Number: R100P68**

REPLACEMENT PART NUMBERS

- OLD R100P09:
 - **NEW Part Number: R100P26** BLENDER RAIL / POLE MOUNT BRACKET
- OLD R100P68:
 - **NEW Part Number: R100P44** HUMIDIFIER RAIL / POLE MOUNT BRACKET

We have made improvements to our Blender and Humidifier brackets to now work on both a pole (15mm-45mm), and a standard hospital rail (10mm x 25mm). We have also made improvements to the knob by using a less bulky, pin-style connection point, which provides better usability when adapting to different diameter poles. **For customers who previously ordered R100P09 or R100P68, please order the above referenced part numbers for all future orders.**



Attached is a datasheet with images of the new mounts.

Kind Regards,

Maxtec Marketing

marketing@maxtec.com


Blender & Humidifier Pole/Rail Mounts

SENSING

ANALYSIS

DELIVERY

*"A single bracket that mounts
to both poles and rails"*



POLE/RAIL MOUNT BRACKETS OVERVIEW

These versatile mounting brackets make attaching blenders and humidifiers to a pole or rail in your hospital more simple and convenient. The design enables a single part to be used on a pole 15mm - 45 mm (.6in to 1.75in), or a standard rail mount (10mm x 25mm). There is a version available for blenders, and a version for humidifiers.

PART NUMBERS

Blender Pole/Rail MountR100P26
Humidifier Pole/Rail Mount . . .R100P44

ML# 0005 REV. A